

# Specification:

#### \* Part

#### 1. Actuator Mechanism

Actuator can be adjust by control that are easily reahable by patient.

## 2. Head & foot panel

Head & foot board made of HDPE, detachable & easy to clean. Equipped with bumper (with humanization design for anti collision effect), located at each corner

#### 3 Wheels

Four wheels 5" double wheel, equipped with brakes on two of them

#### 4. Side guard

Made of aluminium alloy, with collapsible system and can be foled

#### 5. Infuse stand

Made of stainless steel pipe, adjustable height with screw lock. Two hanger (hook).

#### 6. Mattress

Dimension 1900 (L) x 860 (W) x 100 (H) mm





#### \* Position

Back raise lifting angle : 0 - 70° (Approx)

Knee raise lifting angle : 0 - 30° (Approx)

Hi-Low : 490 - 730 mm (Approx)

## \* Dimension (Approx)

Overall 2140 (L) x 970 (W) x 970 (H) mm

#### \* Attention!



Maximum Load 200 Kg

## \* Direction

Push the triangle with upwards point to rise the back, knee and hi-low movement of bed and push the triangle downward point to lower.

#### \* Kemenkes RI

AKD 20902022158

## \* Production Code HBE-LW.3A

#### \* Scan here to visit webpage



